

# FissionGraphed v1.0.8

## User Guide

Tymoteusz Mastalerz\*

2025

### Contents

<b>1</b>	<b>Overview</b>	<b>2</b>
<b>2</b>	<b>Installation and Setup</b>	<b>2</b>
<b>3</b>	<b>Deinstallation process</b>	<b>2</b>
<b>4</b>	<b>Interface Overview</b>	<b>3</b>
<b>5</b>	<b>Getting Started</b>	<b>3</b>
<b>6</b>	<b>Data Visualisation Features</b>	<b>4</b>
<b>7</b>	<b>Advanced Functions</b>	<b>4</b>
<b>8</b>	<b>Menu Reference</b>	<b>4</b>
<b>9</b>	<b>Troubleshooting</b>	<b>5</b>
<b>10</b>	<b>Technical Specifications</b>	<b>5</b>

---

\*[kamilomat12@gmail.com](mailto:kamilomat12@gmail.com)

# 1 Overview

**FissionGraphed** is a scientific visualisation tool that produces interactive scatter plots for nuclear fission data. Primary goals:

- display relationships between fragment and compound–nucleus observables,
- allow quick filtering by publication or atomic number,
- export graphics in PNG format.

## 2 Installation and Setup

### System Requirements

- Operating System: Linux, Debian distribution derivatives if using .deb package variant, any distribution of Linux when using .AppImage variant.
- Python or any other coding language/library is not required.

### **.deb package variant (Debian derivative distributions)**

#### Installation

In order to install the program run command: "sudo dpkg -i ./[package].deb" in terminal while in the same directory as the .deb package.

Alternatively you can run command: "sudo apt install ./[package].deb" in terminal while in the same directory as the .deb package.

#### Running

After installation the application can be found by searching "FissionGraphed" on the App Grid.

### **.AppImage variant (non-Debian distributions)**

#### Installation

In order to use the program you just need to download the .AppImage file.

#### Running

After downloading, go to the file's directory and run it, note, that it won't be present on the App Grid.

## 3 Deinstallation process


### **.deb package variant (Debian derivative distributions)**

In order to uninstall the program you need to run: "sudo apt remove fissiongraphed" in terminal, then go to /opt and remove "fissiongraphed" folder with all it's contents.

### **.AppImage variant (non-Debian distributions)**

In order to uninstall the program you need to delete the .AppImage file.

## 4 Interface Overview



**Plot Area** Four synchronised scatter plots ( $N_L$ ,  $N_H$ ,  $Z_L$ ,  $Z_H$ ).

**Control Panel** Located left; contains publication selector, variant options, axis choice, colour settings and special filters.

**Menu Bar** File, Settings, Filter, View and Help menus for additional actions.

## 5 Getting Started

1. **Select Publications:** Use the multi-select listbox to choose one or more studies. Hold **Ctrl** for additive or subtractive selection.
2. **Choose Variant:**
  - *Inner mode* (1) — to only show inner mode (always present).
  - *Outer mode* (2) — to only show outer mode of fission graph, (shows inner mode for instances with only one mode).
  - *Both modes* — shows both fit modes, inner mode points have solid outline, outer mode points have dashed outlines.
3. **Set X-Axis:** Pick one of six variables; all four subplots update simultaneously.
4. **Customise Appearance:**
  - *View* → *Color by  $Z_{CN}$*  for per-element colouring.
  - *View* → *Filled Markers* for filled or hollow symbols.
  - *Special Chatillon/Bockstiegel* frames for unique styling rules.
5. **Export Plot:** *File* → *Save plot...* to save a PNG (300 dpi). Default filename embeds a timestamp.

## 6 Data Visualisation Features

### Targets

$N_L / N_H$	Neutron numbers of light / heavy fission fragments.
$Z_L / Z_H$	Proton numbers of light / heavy fission fragments.

### Axis Variables

$Z_{CN}$ ,  $N_{CN}/Z_{CN}$ ,  $Z_{CN}^2/A_{CN}$ ,  $N_{CN}-2Z_{CN}$ ,  $N_{CN}$ ,  $A_{CN}$ .

### Colour Schemes

- **Publication-based** (default) — 11 colours  $\times$  11 marker shapes.
- **$Z_{CN}$ -based** — 33 distinct colours; legend displays only present  $Z_{CN}$  values.

### Legends

Automatically rebuild when filters or colouring modes change. Located right of the plots; detached from figure layout.

## 7 Advanced Functions

### Filter by $Z_{CN}$

1. Open *Filter*  $\rightarrow$  *Filter by  $Z_{CN}$* .
2. Select one or more atomic numbers (“All” resets the filter).
3. Click **OK**; plots refresh instantly.

### Set Axis Ranges

1. Choose *Settings*  $\rightarrow$  *Set Axis Ranges*.
2. For each subplot enter **xmin**, **xmax**, **ymin**, **ymax**. Leave blank for auto-scale.
3. Press **Apply**.

### Special Highlighting

<b>Chatillon:Tor</b>	When enabled, Tor nucleus points from Chatillon’s paper appear black and enlarged.
<b>Bockstiegel</b>	Optional gray colouring for all data points from Bockstiegel’s paper.

## 8 Menu Reference

<b>File</b>	<i>Save plot...</i> — export current figure to PNG.
<b>Settings</b>	<i>Set Axis Ranges...</i> — manual axis limits.
<b>Filter</b>	<i>Filter by <math>Z_{CN}</math>...</i> — atomic-number filter.
<b>View</b>	<i>Color by <math>Z_{CN}</math></i> (toggle), <i>Filled Markers</i> (toggle).
<b>Help</b>	<i>Quick Help</i> (popup), <i>User Guide</i> (opens this PDF), <i>About</i> .

## 9 Troubleshooting

### Application Does Not Launch

- Verify `FissionData.dat` is present and readable in `/opt/fissiongraphed/_internal`.
- If problem still persists, best solution is to reinstall the program.

### No Data Visible

- Ensure at least one publication is selected.
- Clear overly restrictive  $Z_{CN}$  filters.

### Export Fails

- Check directory write permissions.
- Confirm sufficient disk space.

## 10 Technical Specifications

### Plot Resolution

PNG, 300 dpi, tight bounding box (legend included).

### Data Format

`.dat` files must be Tab-separated; required columns: `Publication`, `Nucleus`, `Z_CN`, target variables plus pre-computed ratios (see header row of `FissionData.dat` file).

*For feature requests or bug reports, contact the author via email or visit the [project repository](#).*